

## Static friction, sliding friction and rolling friction

### Objects of the experiment

- Investigating the static and sliding friction as a function of the area, the weight and the material
- Comparison of static and sliding friction as a function of the weight and determining the coefficient of friction
- Comparison of rolling and sliding friction as a function of the weight and determining the coefficient of friction

### Principles

In discussing friction between solid bodies, we distinguish between static friction, sliding friction and rolling friction. Static friction force is the minimum force required to set a body at rest on a solid base in motion. Analogously, dynamic (kinetic), or sliding friction force is the force required to maintain a uniform motion of the body. Rolling friction force is the force which maintains the uniform motion of a body which rolls on another body.

To begin, this experiment verifies that the static friction force  $F_H$  and the sliding friction force  $F_G$  are independent of the size of the contact surface and proportional to the resting force  $G$  on the base surface of the friction block. Thus we can say:

$$F_H = \mu_H \cdot G \quad (I)$$

and

$$F_G = \mu_G \cdot G \quad (II).$$

The coefficients of friction  $\mu_H$  and  $\mu_G$  depend on the material of the contact surfaces. As the static friction force is always greater than the sliding friction force, we can always say

$$\mu_H > \mu_G \quad (III).$$

To distinguish between sliding and rolling friction, the friction block is placed on top of multiple stand rods laid parallel to each other. The rolling friction force

$$F_R = \mu_R \cdot G \quad (IV)$$

is measured as the force which maintains the friction block in a uniform motion on the rolling rods. The sliding friction force  $F_G$  is measured once more for comparison, whereby this time the friction block is pulled over the stand rods arranged as a fixed base. The experiment confirms the relationship

$$\mu_G > \mu_R \quad (V).$$

### Setup


- Prepare clean, dry and smooth experiment surfaces (e.g. laboratory bench) for the friction experiments.
- If the resulting frictional forces are too slight, use different base surfaces.

### Carrying out the experiment

- Using the dynamometer, determine the weight (force of gravity)  $G_1$  of the large wooden block and  $G_2$  for the small block.

#### a) Static and sliding friction as a function of the area, the weight and the material

- Place the small block on the experiment surface with the plastic side down.
- Using the dynamometer, measure the maximum horizontal pulling force at which the body remains stationary on the experiment surface as the static friction force  $F_H$ .
- Measure the horizontal pulling force which maintains a uniform motion on the experiment surface as the sliding friction force  $F_G$ .
- Place the wooden block on the base surface with the wide wooden side and then the narrow wooden side down and repeat your measurements for  $F_H$  and  $F_G$ .
- Repeat the measurements with the large block for friction experiments.
- Repeat the measurement on other surfaces as desired.



**Apparatus**

1 Pair of wooden blocks for friction experiments . . . . .	342 10
1 Set 7 weights, 0.1 – 2 kg, with hook . . . . .	315 36
1 Dynamometer, 10.0 N . . . . .	314 47
6 Stand rods 10 cm . . . . .	300 40




Fig. 1 Measuring the static friction force  $F_H$  (top) and the sliding friction force  $F_G$  (bottom).




Fig. 2 Measuring the static friction force  $F_H$  (top) and the sliding friction force  $F_G$  (bottom) as a function of the force of gravity  $G$ .




Fig. 3 Measuring the rolling friction force  $F_R$  (top) and the sliding friction force  $F_G$  (bottom) as a function of the force of gravity  $G$ .

**b) Static and sliding friction force as a function of the force of gravity**

- Place the large block on the experiment surface with the plastic side down and measure the static and sliding friction force.
- Increase the weight of the block by adding in turn the weights 0.1 kg, 0.2 kg, 0.5 kg and 1.0 kg and repeat the measurements.
- Carry out the same measurements for the wooden side of the block as well.

**c) Rolling and sliding friction force as a function of the force of gravity**

- Lay the stand rods next to each other and place the large block on the rods with the plastic side down.
- Measure the horizontal pulling force which maintains a uniform motion on the rolling rods as the rolling friction force  $F_R$ .
- Increase the weight of the block by adding in turn the weights 0.1 kg, 0.2 kg, 0.5 kg and 1.0 kg and repeat the measurements.
- Align the block parallel to the rod axes and measure the sliding friction force.

**Measuring example**

$G_1 = 1.5 \text{ N}$

$G_2 = 3.0 \text{ N}$

**a) Static and sliding friction as a function of the area, the weight and the material**

Experiment surface: plastic-coated benchtop

Tab. 1: Static friction force  $F_H$  and sliding friction force  $F_G$  as a function of the force of the area, the weight and the material

$\frac{G}{N}$	Material	$\frac{A}{\text{cm}^2}$	$\frac{F_H}{N}$	$\frac{F_G}{N}$
1.5	Plastic	$12 \times 6$	0.8	0.6
1.5	Wood	$12 \times 6$	0.3	0.3
1.5	Wood	$12 \times 3$	0.3	0.3
3.0	Plastic	$12 \times 6$	1.6	1.1
3.0	Wood	$12 \times 6$	0.5	0.5

**b) Comparison of static and sliding friction force**

Experiment surface: plastic-coated benchtop

Table 2: Static friction force  $F_H$  and sliding friction force  $F_G$  as a function of the force of gravity  $G$

$\frac{G}{N}$	Plastic side		Wooden side	
	$\frac{F_H}{N}$	$\frac{F_G}{N}$	$\frac{F_H}{N}$	$\frac{F_G}{N}$
3	1.6	1.1	0.5	0.4
4	2.2	2.0	0.8	0.6
5	3.1	2.8	0.9	0.8
8	5.0	4.6	1.9	1.3
13	8.3	8.0	3.0	2.0

c) Comparison of sliding and rolling friction

Table 3: Sliding friction force  $F_G$  and rolling friction force  $F_R$  as a function of the force of gravity  $G$

$\frac{G}{N}$	$\frac{F_G}{N}$	$\frac{F_R}{N}$
3	3.0	0.1
4	4.3	0.2
5	5.2	0.2
8	9.0	0.3
11		0.4
13		0.5
18		0.6
23		0.7

Evaluation and results

a) Static and sliding friction as a function of the area, the weight and the material

As the measuring results in Table 1 show, both the static friction force and the sliding friction force depend on the material properties of the friction surfaces and on the weight (force of gravity) of the blocks. However, the friction forces are independent of the size of the friction area.

b) Comparison of static and sliding friction force




Fig. 4 Static friction force  $F_H$  and sliding friction force  $F_G$  as a function of the force of gravity  $G$

- Closed squares: static friction for plastic
- Open squares: sliding friction for plastic
- Closed circles: static friction for wood
- Open circles: sliding friction for wood
- Experiment surface: plastic-coated benchtop

Fig. 4 shows the results of the measurements. The slope of the line through the origin is equivalent to the coefficients of friction  $\mu_H$  and  $\mu_G$  calculated using (I) and (II) (see Table 2).

Table 4: Static friction coefficient  $\mu_H$  and sliding friction coefficient  $\mu_G$  for friction on plastic

Material	$\mu_H$	$\mu_G$
Plastic	0.63	0.59
Wood	0.22	0.15

c) Comparison of sliding and rolling friction




Fig. 5 Sliding friction force  $F_G$  (circles) and rolling friction force  $F_R$  (triangles) as a function of the force of gravity  $G$

Fig. 5 shows the results of the measurements. The slope of the line through the origin is equivalent to the coefficients of friction  $\mu_H$  and  $\mu_G$  calculated using (I) and (II) (see Table 3).

Table 5: Rolling friction coefficient  $\mu_R$  and sliding friction coefficient  $\mu_G$

Material:	$\mu_G$	$\mu_R$
Plastic	1.09	0.03