


Acoustic beats



can also be carried out with [Pocket-CASSY](#)

■ Load example

Experiment description

This experiment records a beat caused by two tuning forks having slightly different frequencies. The individual frequencies f_1 and f_2 , the new oscillation frequency f_n and the beat frequency f_s are determined and compared with the theoretical values

$$f_n = \frac{1}{2} (f_1 + f_2) \quad \text{and} \quad f_s = |f_1 - f_2|.$$

Equipment list

1	Sensor-CASSY	524 010
1	CASSY Lab	524 200
1	Pair of resonance tuning forks	414 72
1	Multi-purpose microphone with Base	586 26 300 11
	or	
1	Microphone S	524 059
1	PC with Windows 98/2000/XP/Vista	

Experiment setup (see drawing)

Position the two tuning forks and the multipurpose microphone and connect the latter to input A of Sensor-CASSY (don't forget to set the microphone to "Signal" mode and switch it on). "Mistune" one of the tuning forks slightly by attaching an additional weight.

Carrying out the experiment

■ Load settings

- Strike the first tuning fork and start the measurement with **F9**.
- Optimize the signal strength using the control on the microphone.
- Determine the frequency f_1 (e. g. by marking the **Standard** display with [vertical lines](#) or finding the [peak center](#) in the **Frequency Spectrum**).



- Delete the measurement with **F4**.
- Strike the second tuning fork and start the measurement with **F9**.
- Find frequency f_2 .
- Delete the measurement with **F4**.
- Strike both tuning forks equally hard and start the measurement with **F9**.

Evaluation

When the amplitudes of the two tuning forks are equal, the nodes and antinodes in the beat become clearly distinguished. The beat frequency f_s is determined from the interval T_s between two nodes as $f_s = 1/T_s$.

To ensure satisfactory accuracy in determining the new oscillation period T_n , you should average the measurements over 10 periods before determining the new oscillation frequency as $f_n = 1/T_n$. To determine the time differences you can e.g. insert [vertical lines](#) in the diagram or measure the [difference](#) directly.

In this example we obtain $f_1 = 425$ Hz, $f_2 = 440$ Hz, $f_n = 433$ Hz, $f_s = 14.5$ Hz, which closely confirms the theory $f_n = \frac{1}{2}(f_1 + f_2) = 432.5$ Hz and $f_s = |f_1 - f_2| = 15$ Hz.

Click on the **Frequency Spectrum** to read out the two frequencies f_1 and f_2 of the tuning forks and their amplitudes. The easiest way to determine the frequency is to find the [peak centers](#).